

## Aspire ComfiMotion Activ Care Bed



The Activ Care bed combines the best of in-home care functionality, positioning and low height safety with luxurious and contemporary style for enhanced comfort, care and wellbeing.

- **Ultra Low Height:** Minimum platform height of 195 mm reduces falls risk, improving user safety
- **Customisable Safe Entry/Exit Height:** Programmable single button press ensures perfect bed height to maximise independence
- **Under Bed Clearance:** provides clear access for patient lifters and overbed tables
- **Battery Backup:** Provides confidence and maintains independence during power outages
- **Night Assist:** Luminous underbed halo for safe entry and exit
- **Illuminated Handset:** Large button interface is intuitive and highly visible at all times
- **Massage:** Relaxing vibration modes for personalised neck, shoulder and back relief
- **Pull-apart Design:** Split frame is more accessible in homes with restricted access such as tight staircases or without a commercial lift.
- **Hi-Lo Functionality:** Assists with manual handling and helps to mitigate fall risks
- **Head Tilt:** Relaxes neck muscles, providing optimal support whilst reading or watching TV
- **Back Raise:** Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring and sleep apnoea
- **Auto-Regression in the Backrest:** Stops the user sliding down the bed during movement while improving comfort and maximising pressure redistribution
- **Clinical Positioning:** One touch raises back and legs to alleviate cardiac



Contract approved

### Support @ Home

**Function:** Beds and bed equipment

**Category:** Prescribed

Click [here](#) for more information

### Compliance

AS/NZS 60601-2-52:2009

exertion, improving resting comfort

- **Leg Raise:** Promotes circulation, relieving tired, swollen/ fatigued or painful limbs
- **Memory Function:** Allows the user to set and return to their preferred positions at any time with one click of a button
- **Zero Gravity:** Helps to relax the circulatory system and muscles of the lower back
- **Seated Position:** Assists with bed entry and exit and upright dining
- **Trendelenburg (and Reverse):** Can assist care givers when repositioning a user\*
- **Lockable Castors:** Eight heavy duty dual castors
- Headboard and accessories sold separately

\*Carers should seek advice from a trained Allied Health professional for correct manual handling techniques.

## Warranty Information

- 10 Years Limited

| Code      | Display Name                                                                      | Total Length | Total Width | Total Height | Safe Working Load | Colour          |
|-----------|-----------------------------------------------------------------------------------|--------------|-------------|--------------|-------------------|-----------------|
| BEB046745 | Aspire ComfiMotion Activ Care Bed - Portobello Grey - Long Single (Headboard Sep) | 2030mm       | 920mm       | 200 - 775mm  | 250kg             | Portobello Grey |
| BEB046755 | Aspire ComfiMotion Activ Care Bed - Portobello Grey - King Single (Headboard Sep) | 2030mm       | 1070mm      | 200 - 775mm  | 250kg             | Portobello Grey |
| BEB046765 | Aspire ComfiMotion Activ Care Bed - Portobello Grey - Long Double (Headboard      | 2030mm       | 1370mm      | 200 - 675mm  | 360kg             | Portobello Grey |

| Code      | Display Name                                                                | Total Length | Total Width | Total Height | Safe Working Load | Colour          |
|-----------|-----------------------------------------------------------------------------|--------------|-------------|--------------|-------------------|-----------------|
|           | Sep)                                                                        |              |             |              |                   |                 |
| BEB046775 | Aspire ComfiMotion Activ Care Bed - Portobello Grey - Queen (Headboard Sep) | 2030mm       | 1520mm      | 200 - 675mm  | 360kg             | Portobello Grey |