

Aspire Deluxe Bedside Commode



- Allows for toileting in the bedroom for users unable to travel between the bedroom and bathroom
- Available in Champagne or Slate Vinyl
- Easy Clean Vinyl Surfaces
- Removable Bowl with Lid & Carry Handle
- Height Adjustable
- Raised Padded Armrests
- Fold Away, Padded Seat
- Non-Slip Tips
- Toilet Seat to Floor Height 435 - 585 mm

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Toileting

Category: Prescribed

Click [here](#) for more information

Warranty Information

- 1 Year

NDIS

Item Code:

05_091203048_0103_1_2

Toilet Aids – Composite:

05_414_0103_1_2

Click [here](#) for more information

Compliance

AS/NZS 3973:2009



Code	Display Name	Seat Depth	Seat Width	Seat To Floor Height	Safe Working Load	Unit Weight
BEC045115	Aspire Deluxe Bedside Commode - Champagne	450 mm	460 mm	410 - 565 mm	160 kg	12 kg
BEC045116	Aspire Deluxe Bedside Commode - Slate	450 mm	460 mm	410 - 565 mm	160 kg	12 kg

