

## Aspire Autoclavable Bed Pan



Our Autoclavable Bed Pans are made from premium polypropylene. They will last for many years in any healthcare environment subject to the rigours of cleaning, sanitation and infection control requirements.

- Heavy duty and reusable
- Designed for use in healthcare environments such as surgeries, acute care, general ward, medical centres, and aged care facilities
- Autoclavable polypropylene that can be cleaned using autoclaves and sanitisers
- Pastel green colour for easy identification

### Warranty Information

- 1 Year



✓ Contracted DVA product

✓ Contract approved

### Support @ Home

**Function:** Collecting urine and faeces

**Category:** Prescribed

Click [here](#) for more information

### NDIS

**Item Code:**

05\_098800044\_0103\_1\_2

**Toilet Aids – Composite:**

05\_414\_0103\_1\_2

Click [here](#) for more information

### Compliance

AS/NZS 3973:2009

| Code      | Display Name | Total Length | Capacity |
|-----------|--------------|--------------|----------|
| BTT151800 | Bed Pan      | 360 mm       | 1000 ml  |

