

## Lateral Chest Support Assembly

- A great range of accessories to fit the Aspire Manual Wheelchair range to improve functionality and user experience
- Click [here](#) to ensure this accessory fits your specific wheelchair model



### Warranty Information

- 1 Year



| Code      | Display Name                                                  |
|-----------|---------------------------------------------------------------|
| WCA848430 | Lateral Chest Support Assembly - Left                         |
| WCA848450 | Rehab RS/RX Wheelchair Lateral Chest Support Assembly - Left  |
| WCA848440 | Lateral Chest Support Assembly - Right                        |
| WCA848460 | Rehab RS/RX Wheelchair Lateral Chest Support Assembly - Right |

✓ Contract approved

### Support @ Home

**Function:** Accessories for wheelchair seating

**Category:** Prescribed

Click [here](#) for more information

### Compliance

AS3696-1:2008